

Laboratory Freezers

Model	ISU27	ISU37	ISU57	ISU97
SKU number	444448572	444449327	444448522	444448523
Product dimensions (in mm)				
Width	500	495	558	595
Depth	510	580	600	600
Height	562	830	863	1550
Packed dimensions (in mm)				
Width	520	560	600	625
Depth	595	655	610	700
Height	590	850	900	1610
Product weight (in Kg)	17	33	38	75
Packed weight (in Kg)	20	36.5	39.5	77.5
Net capacity (in litres)	38	87	94	202
Gross capacity (in litres)	50	102	116	230
Temperature range	-18°C to -20°C	-18°C to -20°C	-18°C to -25°C	-18°C to -25°C
Colour	White	White	White	White
Door type	Hinged, Reversible	Hinged, reversible	Hinged, reversible	Hinged, reversible
Door style	Solid	Solid	Solid	Solid
Replaceable door seal	-	-	-	Yes
Open door sensor	-	-	-	-
Lock	Yes	Yes	Yes	Yes
Adjustable feet	Fixed	Yes	Yes	Yes
	Fixed	Fixed	Fixed	Roller Castors
Temperature display in the panel	Yes	Yes	Yes	Yes
Temperature alarm	Yes	Yes	Yes	Yes
Shelves	Yes (1)	Yes (2)	-	-
Wire drawer	-	Yes (1)	Yes (4)	Yes (6)
Interior light	-	-	-	-
Freezer Defrost	Manual	Manual	Manual	Manual
Fan	-	-	-	-
Control	Electronic	Electronic	Electronic	Electronic
Climate class	N	N	N	N
Coolant	R600a	R134a	R600a	R600a
AUX Output	-	-	-	-
UK 3-pin plug	Yes	Yes	Yes	Yes
Voltage supply	230V / 50Hz	230V / 50Hz	230V / 50Hz	230V / 50Hz