

# Laboratory Fridge-Freezers

Model	ISR27	ISR57	IST47
SKU number	444448571	444448518	444448520
Product dimensions (in mm)			
Width	502	558	560
Depth	540	600	600
Height	660	863	1572
Packed dimensions (in mm)			
Width	520	600	600
Depth	595	610	610
Height	675	900	1620
Product weight (in Kg)	25	36	68
Packed weight (in Kg)	27	37.5	70
Net capacity (in litres)	75	110	152
	Fridge		
	Freezer		
Gross capacity (in litres)	83.5 (77/6.5)	130 (112/18)	246 (157/89)
Temperature range	+2°C to +10°C / -6°C to -9°C	+2°C to +10°C / min -18°C	+2°C to +10°C / 18°C to -25°C
Colour	White	White	White
Door type	Hinged, reversible	Hinged, reversible	Hinged, reversible
Door style	Solid	Solid	Solid
Replaceable door seal	-	-	-
Open door sensor	-	-	-
Lock	Yes	Yes	Yes (2)
Adjustable feet	Fixed	Yes	Yes
	Front		
	Back		
Temperature display in the panel	Fixed	Fixed	Roller Castors
Temperature alarm	Yes	Yes	Yes (Fridge)
Shelves	Yes (2)	Yes (2)	Yes (3) fridge
			Yes (1) fridge /
Wire drawer	Yes (1)	-	Yes (3) freezer
Moulded drawer	-	Yes (1)	-
Adjustable Door Bins	-	Yes (3)	-
Interior light	-	-	-
Fridge defrost	Automatic	Automatic	Automatic
Freezer Defrost	Manual	Manual	Manual
Fan			
	-	-	-
Control	Electronic	Electronic	Electronic
Climate class	N	N	N
Coolant	R134a	R600a	R600a
AUX Output	-	-	-
UK 3-pin plug	Yes	Yes	Yes
Voltage supply	230V / 50Hz	230V / 50Hz	230V / 50Hz