

# Laboratory Fridges

Model	LR207	LR307	LR507	LR907	LR1607
SKU number	444448568	444449326	444448514	444448516	444448569
Product dimensions (in mm)					
Width	502	495	558	595	685
Depth	540	580	600	600	910
Height	660	830	863	1550	1688
Packed dimensions (in mm)					
Width	520	560	600	625	770
Depth	595	655	610	700	950
Height	675	850	900	1610	1810
Product weight (in Kg)	27.5	31.5	34	81.5	130
Packed weight (in Kg)	30.5	35	35.5	84	140
Net capacity (in litres)	67	103	145	285	441
Gross capacity (in litres)	82	107	153	286	475
Temperature range	+2°C to +10°C	+2°C to +10°C	+2°C to +10°C	+2°C to +10°C	+2°C to +10°C
Colour	White	White	White	White	White
Door type	Hinged, reversible	Hinged, reversible	Hinged, reversible	Hinged, reversible	Hinged R/H only
Door style	Solid	Solid	Solid	Solid	Solid
Replaceable door seal	-	-	-	Yes	-
Open door sensor	-	-	-	-	-
Lock	Yes	Yes	Yes	Yes	Yes
Adjustable feet	Fixed	Yes	Yes	Yes	Lockable Roller Castors
		Fixed	Fixed	Fixed	Roller Castors
Temperature display in the panel	Yes	Yes	Yes	Yes	Yes
Temperature alarm	Yes	Yes	Yes	Yes	Yes
Shelves	Yes (2)	Yes (2)	Yes (2)	Yes (5)	Yes (5)
Wire drawer	Yes (1)	Yes (1)	Yes (1)	Yes (1)	-
Interior light	-	-	-	-	-
Fridge defrost	Automatic	Automatic	Automatic	Automatic	Automatic
Fan	-	-	-	-	-
Control	Electronic	Electronic	Electronic	Electronic	Electronic
Climate class	N	N	N	N	N
Coolant	R134a	R134a	R600a	R600a	R134a
AUX Output	-	-	-	-	-
UK 3-pin plug	Yes	Yes	Yes	Yes	Yes
Voltage supply	230V / 50Hz	230V / 50Hz	230V / 50Hz	230V / 50Hz	230V / 50Hz