

# Pharmacy Glass Door

Model	PG207	PG507	PG907	PG1607
SKU number	444448565	444448509	444448511	444448566
<b>Product dimensions (in mm)</b>				
Width	502	558	595	685
Depth	540	600	600	910
Height	660	863	1550	1668
<b>Packed dimensions (in mm)</b>				
Width	520	600	625	770
Depth	595	610	700	950
Height	675	900	1610	1810
Product weight (in Kg)	25	52	107	130
Packed weight (in Kg)	27	54	110	140
Net capacity (in litres)	67	137	275	393
Gross capacity (in litres)	82	153	286	444
Temperature range	+2°C - +8°C	+2°C - +8°C	+2°C - +8°C	+2°C - +8°C
Colour	White	White	White	White
Door type	Hinged, reversible	Hinged, reversible	Hinged, reversible	Hinged R/H only
Door style	Glass	Glass	Glass	Glass
Replaceable door seal	Yes	Yes	Yes	Yes
Open door sensor	Yes	Yes	Yes	Yes
Lock	Yes	Yes	Yes	Yes
Adjustable feet	Front	Fixed	Yes	Yes
	Back	Fixed	Fixed	Roller Castors
Temperature display in the panel	Yes	Yes	Yes	Yes
Temperature alarm	Yes	Yes	Yes	Yes
Shelves	Yes (2)	Yes (2)	Yes (5)	Yes (5)
Wire drawer	Yes (1)	Yes (1)	Yes (1)	-
Interior light	Yes (6W)	Yes (15W)	Yes (15W)	Yes (12W)
Fridge defrost	Automatic	Automatic	Automatic	Automatic
Fan	Yes	Yes	Yes	Yes
Control	Electronic	Electronic	Electronic	Electronic
Climate class	N	N	N	N
Coolant	R134a	R600a	R600a	R134a
AUX Output	-	Yes	Yes	Yes
UK 3-pin plug	Yes	Yes	Yes	Yes
Voltage supply	230V / 50Hz	230V / 50Hz	230V / 50Hz	230V / 50 Hz