



Under-counter Laboratory Spark Free Freezer 98L

GDPA code: **LSFSF98UK**

SKU code: **444410535**

- ✓ -18°C to -25°C temperature range
- ✓ Spark free interior
- ✓ External digital controller and temperature display
- ✓ Fitted lock and 2 keys
- ✓ Audible and visual temperature alarms



The Experts in
Medical Refrigeration

www.lec-medical.co.uk
0344 815 3742


1 of 3



Specification: Under-counter Laboratory Spark Free Freezer 98L

Gross Capacity (litres)	98
Net Capacity (litres)	84
Unit External Width (mm)	595
Unit External Depth (mm)	642
Unit External Height (mm)	850
Product Weight (kg)	37.25
Packed Width (mm)	640
Packed Depth (mm)	660
Packed Height (mm)	910
Packed Weight (kg)	41.75
Country of Origin	UK
Internal Dimensions (mm) (WxDxH)	420 x 430 x 640
Temperature Range	-18°C to -25°C (set point tolerance + or - 2)
Colour	White/Green
External Finish	Painted Steel
Internal Finish	HIPS*
Type of plastic	SRL800 (HIPS/ABS mix)
Door Finish	Painted Steel / Solid
Door Type	Solid
Door Style	Hinged, Reversible
Door Open Depth (mm)	1130
Replacement Door Seal	No
Open Door Sensor	Yes
Lock	Yes, 2 Keys
Adjustable Feet (number & position)	4 - All adjustable
Porthole	Yes
Porthole Size	16 mm
Porthole Location	Rear LH Corner
Temperature Display on Panel	Yes
Temperature Alarm	Yes
Remote Alarm Contacts	Yes
Audible Alarms	Yes
Visual Alarms	Yes
Number of Shelves	3
Number of Shelf Positions	3
Shelf Size (mm) (WxDxH)	495 x 325
Shelf Material	Glass / Plastic
Number of Baskets	0
Basket Material	N/A
Basket Size (mm) (WxDxH)	485 x 180 x 171
Interior Light	No
Light Type	N/A
Light Activation	N/A
Fan	No
Spark Free	Yes
Auto or Manual Defrost	None
Controller Type	Digital, External
Number of Probes	1
Probe Type	Air
Climate Class	SN
Refrigerant Type	R600a

Specification: Under-counter Laboratory Spark Free Freezer 98L



CFC Free	Yes
HFC Free	Yes
HCFC Free	Yes
Ammonia Free	Yes
UK 3 Pin Plug	Yes
Other Plug Types Available	Yes
Voltage Supply	230v
Insulation Thickness (mm)	85
Insulation Type	High Density Polyurethane Foam
Blowing Agent	Cyclopentane
Anti Microbial Handle	Yes
Optional Drawer System Available	No
Wall Brackets	No
Digital Lock	No