



Solid Door

Glass Door

Under-counter Control Plus Solid/Glass Door Refrigerator 158L

GDPA code: PPSR158UK (solid), PPGR158UK (glass)

SKU code: 444410526 (solid), 444410527 (glass)


- ✓ 2°C to 8°C temperature range
- ✓ Air and core temperature display
- ✓ Represents your vaccine vial
- ✓ External digital controller and temperature display
- ✓ Fitted lock and 2 keys
- ✓ Audible and visual temperature alarms
- ✓ 5 year warranty as standard



Specification: Under-counter Control Plus Solid/Glass Door Refrigerator 158L

Gross Capacity (litres)	158
Net Capacity (litres)	154
Unit External Width (mm)	595
Unit External Depth (mm)	642
Unit External Height (mm)	850
Product Weight (kg)	33.25 (solid), 36.25 (glass)
Packed Width (mm)	640
Packed Depth (mm)	660
Packed Height (mm)	910
Packed Weight (kg)	37.82 (solid), 39.75 (glass)
Country of Origin	UK
Internal Dimensions (mm) (WxDxH)	500 x 450 x 725. Depth of Step 300 (Door Closed)
Temperature Range	+2°C to +8°C (tolerance + or - 2)
Colour	White/Green
External Finish	Painted Steel
Internal Finish	HIPS*
Type of plastic	SRL800 (HIPS/ABS mix)
Door Finish	Painted Steel / Solid
Door Type	Solid (solid), Glass (glass)
Door Style	Hinged, Reversible
Door Open Depth (mm)	1130
Replacement Door Seal	No
Open Door Sensor	Yes
Lock	Yes, 2 Keys
Adjustable Feet (number & position)	4 - All adjustable
Porthole	Yes
Porthole Size	16 mm
Porthole Location	Rear LH Corner
Temperature Display on Panel	Yes
Temperature Alarm	Yes
Remote Alarm Contacts	Yes
Audible Alarms	Yes
Visual Alarms	Yes
Number of Shelves	3
Number of Shelf Positions	6
Shelf Size (mm) (WxDxH)	496 x 405
Shelf Material	Wire
Number of Baskets	1
Basket Material	Wire
Basket Size (mm) (WxDxH)	485 x 255 x 171
Interior Light	Yes
Light Type	LED
Light Activation	Door Switch / Controller
Fan	Yes - Tangential
Spark Free	No
Auto or Manual Defrost	Both (solid), None (glass)
Controller Type	Digital, External
Number of Probes	2
Probe Type	Air & Load (Glycol Solution)
Climate Class	SN
Refrigerant Type	R600a

Specification: Under-counter Control Plus Solid/Glass Door Refrigerator 158L



CFC Free	Yes
HFC Free	Yes
HCFC Free	Yes
Ammonia Free	Yes
UK 3 Pin Plug	Yes
Other Plug Types Available	Yes
Noise Level (dB)	TBC
Energy Consumption (kwh)	TBC
Voltage Supply	230v
Aux Output	No
Insulation Thickness (mm)	55
Insulation Type	High Density Polyurothane Foam
Blowing Agent	Cyclopentane
Anti Microbial Handle	Yes
Optional Drawer System Available	No
Wall Brackets	No
Digital Lock	No